

# Let's talk!

## ...about safety gear for LN2 handling

Protect yourself from cold burns when handling extreme cold objects. Liquid nitrogen is below  $-196^{\circ}\text{C}$  and may cause severe skin damage.

Recommended protection:

- Cryogenic gloves
- Face shield or goggles
- Cryogenic apron
- Long sleeves, cuffless pants
- Closed toe shoes



Model and size	Wrist (30 cm)	Mid arm (37 cm)	Elbow (60 cm)
Small 8	GLO/CMS/KW	GLO/CMS	GLO/CMES
Medium 9	GLO/CMM/KW	GLO/CMM	GLO/CMEM
Large 10	GLO/CML/KW	GLO/CML	GLO/CMEL
Price (vat 0%)	235 €	235 €	275 €



Alcohol and water-resistant cryogenic marker

- Marked surfaces may be exposed to LN2 or  $-80^{\circ}\text{C}$
- Permanent, black marker
- Fine tip
- Ink is water-resistant and alcohol resistant



Handheld printer for lab printing needs. Nylon and self-laminating polyester label is suitable for LN2 application, freezer and hot water

Cat no	Description	Price (vat 0%)
M210-LAB kit EU	Label printer kit packed in a hard case containing additionally Li-Ion battery, powercord, label cartridge M21-750-7425	260 €
M21-375-461-AW	Self-laminating polyester cassette width 9,53 mm	65 €
M21-500-461-AW	Self-laminating polyester cassette width 12,7 mm	75 €
M21-750-461-AW	Self-laminating polyester cassette width 19,05 mm	85 €
M21-375-499	Nylon cassette width 9,53 mm	48 €
M21-500-499	Nylon cassette width 12,7 mm	52 €
M21-750-499	Nylon cassette width 19,05 mm	65 €

*Valid in February 2023 with promo code H223*

📞 020 7410 270

✉ info@bionordika.fi

🏠 shop.bionordika.fi

