

Special Offers

for Selected Corning Equipment



Axygen® Mini Plate Spinner Centrifuge

- Capacity for two standard plates
- Top loading design
- Use before and after thermal cycling to improve PCR yield
- Quickly spin down droplets and condensation
- Accepts skirted and non-skirted and all standard PCR plates

Tarjoushinta 495€
(listahinta 717€)

Cat.no. PLATESPINNER-230EU



Falcon® Pipet Controller

- Designed for use with glass and plastic serological pipets
- Aspirating and dispensing speed is controlled through the finger triggers.
- the batteries enable up to 8 hours of continuous use
- large LCD display clearly indicates battery status, pipetting mode and speed.

Tarjoushinta 175€
(listahinta 341€)

Cat. no. 357469



Corning® LSE™ Low Speed Orbital Shaker

- Ideal for staining and washing of fragile gels
- Gentle motion prevents foaming of liquids
- includes a 30 x 30 cm platform
- Speed can be set from 3 to 60 rpm
- Simple analog control with easy-turn knobs
- 19mm orbit for blotting, gel staining and general mixing
- Safe for cold room and incubator use

Tarjoushinta 850€
(listahinta 1.107€)

Cat. no 6781-FP



Corning® LSE™ Digital Vortexer

- Three versatile platform options (available separately)
- Compact design with easy-access control panel
- Operation of the vortexer can be timed or continuous, while the speed is variable.
- Non-slip rubber mat to keep containers in place
- Safe for cold room or incubator use

Tarjoushinta 840€
(listahinta 1.207€)

Cat. no. 6710



Corning® LSE™ Platform Rocker, Single Platform

- high capacity rocker that is ideal for western blotting and general mixing applications
- Non-slip rubber pads and 5 kg load capacity
- adjustable rocking speed from 2 to 60 rpm and fixed 7 degree tilt
- safe to use in temperature controlled environments up to 70°C
- single platform

Tarjoushinta 650€
(listahinta 882€)

Cat. no. 6704



Corning® LSE™ Platform Rocker, Double Platform

- high capacity rocker that is ideal for western blotting and general mixing applications
- Non-slip rubber pads and 5 kg load capacity
- adjustable rocking speed from 2 to 60 rpm and fixed 7 degree tilt
- safe to use in temperature controlled environments up to 70°C
- double platform

Tarjoushinta 690€
(listahinta 967€)

Cat. no. 6706