

NEW!



BIOTOPICS

Biotop NEWS & PROMOTIONS END-OF-YEAR 2022

Protein Electrophoresis & Western Blot

FastGene Western Blot System

Wet-transfer Western Blot system. Runs two blots simultaneously.

FastGene Western ECL Kit

Western blot detection kit, high sensitivity. Easy switch from other brands.

FastGene Block & Go

Western blocking solution and signal enhancer. 2-in-1 solution: Blocking reagent and signal enhancer.



*Nippon Genetics:
Protein Electrophore-
sis, Western Blotting*



FastGene PAGE Protein System

Complete set, compatible with different pre-cast gels.

FastGene Q-Stain

Ultra-fast protein staining reagent - results in 10 minutes or less. Single-Step procedure – No washing, fixing or destaining.

BlueStar Protein Markers

Prestained three colour protein ladders. 10-12 sharp bands with 3 colours (blue, red and green).

 **NIPPON
Genetics**

ORDERING INFORMATION

		PRICE
PG01	FastGene PAGE Protein System	1480 €
PG08	FastGene Western Blot System	815 €
FG-CH05	FastGene Block & Go	98 €
FG-CH01	ECL Kit. Western blot detection kit (100 ml)	58 €
MWP04	BlueStar Plus Prestained Protein Marker (500 µl)	120 €
FG-QS1	FastGene Q-Stain (1L)	189 €



BIOTOP JL

www.biotop.fi

Prices valid until the end of December 2022. Vat 0%.

Detachment Solutions

The DetachKits are easy to use and contain all components needed for enzymatic cell detachment.



DetachKit

For adherent primary human/animal cells and established cell lines. Also suitable for serum free cell culture.

DetachKit 2

For more sensitive cell types, e.g. Promo-Cell HNEpC.

Suitable for serum free cell culture.

ORDERING INFORMATION

DetachKit

C-41200	3 x 30 ml	87 €
C-41210	3 x 125 ml	244 €
C-41220	3 x 250 ml	409 €

DetachKit 2

C-41202	3 x 30 ml	90 €
C-41212	3 x 125 ml	266 €
C-41222	3 x 250 ml	449 €

MADE IN GERMANY



Low DNA Binding Plates

For Next Generation Sequencing & Other Sensitive Applications

- Maximum nucleic acid sample recovery after incubation
- No chemical additives or coatings used to achieve low binding characteristics
- Better plastic consumables for wider use in more applications
- Several suitable plate formats available (other formats available on request)

ORDERING INFORMATION

FrameStar Low Binding PCR Plates (50 plates)

4ti-LB0770/C FrameStar 96 Semi-Skirted, ABI style

4ti-LB0960/RIG FrameStar 96 Skirted, low profile Plate, extra rigid

4ti-LB0384/RIG FrameStar 384 Skirted, low profile Plate, extra rigid

Low Binding Microplates (50 plates)

4ti-LB0125 96 Deep Well Storage Microplate, for magnetic separators

4ti-LB0109 96 Round Well Microplate,

4ti-LB0125 96 Round Deep Well Storage Microplate, for magnetic separators

MADE IN GREAT BRITAIN



AZENTA
LIFE SCIENCES



4titude has been rebranded and is now part of the
Azenta Life Sciences portfolio

Prices valid until the end of December 2022. Vat 0%.



MACHEREY-NAGEL

Classics Purification Kits

Promotion prices
plus Fazer Chocolate!

Isolation of RNA, plasmid DNA
and genomic DNA, PCR Clean-up
and Gel Extraction,
Purification of His-tag Proteins



MACHEREY-NAGEL
Classics Purification
Kits. Promotion Prices!

MADE IN GERMANY



FOR SYBR® GREEN AND PROBE ASSAYS

SensiFAST™ qPCR Kits

Ultra-sensitive and specific kits for a fast, broad dynamic range and highly reproducible qPCR.

SensiFAST™ qPCR mixes are designed for all qPCR platforms and they are produced, unlike other common qPCR mixes, **under ISO 13485** which makes them suitable to use as an IVD component. SensiFAST™ Probe is designed for TaqMan® probes.

ORDERING INFORMATION **PROMOTION PRICES!**

SensiFAST™ SYBR®, 2000 rxn (4x5 ml) (20x5 ml)

BIO-92020 Hi-ROX Kit **655€ 2975€**

BIO-94020 Lo-ROX Kit **655€ 2975€**

BIO-98020 No-ROX Kit **655€ 2975€**

SensiFAST™ Probe, 2000 rxn

BIO-82020 Hi-ROX Kit **655€ 2975€**

BIO-84020 Lo-ROX Kit **655€ 2975€**

BIO-86020 No-ROX Kit **655€ 2975€**

ALSO AVAILABLE!

**SensiFAST™ Probe
One-Step Kit for fast
RT-qPCR.**

INSTRUMENT COMPATIBILITY chart at www.biotop.fi: TECHNICAL

SENSITIVE & FAST

SensiFAST™ cDNA Synthesis Kit

Starting from 1 pg total RNA. Highly sensitive 1st strand synthesis with unbiased and robust efficiency leading to enhanced reproducibility and data accuracy.

ORDERING INFORMATION

SensiFAST™ cDNA Synthesis Kit

BIO-65053 50 rxn **198€**

BIO-65054 250 rxn **649€**

meridian BIOSCIENCE™
LIFE DISCOVERED. LIFE DIAGNOSED.

Bioline products are produced under ISO 13485 and are suitable for further manufacturing use as an IVD component. **MADE IN GERMANY.**

mRNA transfection reagent

jetMESSENGER®

High efficiency on a wide variety of difficult
to transfect cells. **REQUEST A QUOTE!**

Polyplus+
transfection™

MADE IN FRANCE





Posti Green



Prices valid until the end of December 2022. Vat 0%.

NOW 985 €
Total 1345 €

Mupid ONE PROMOTION!

Simply the best DNA Electrophoresis system!

Mupid ONE + the Blue LED LID + MIDORI Green XTRA

ALL YOU NEED!

Mupid Bundle - PROMOTION!

MU2 Mupid-One Elektrophoresis System **843 €**

MU4 Mupid-One LED Illuminator **398 €**

MG10 Midori Green Xtra (1 ml) **104 €**

MADE IN JAPAN

FastGene® Instruments PROMOTION!



Mini Centrifuge

Available in four different colours: blue, green, red and pink (NG002)

Vortexer Mini
Continuously adjustable vortexer (VX2)



Plate Centrifuge

For 96- and 384-well plates (NG040)



Mini Dry Bath Advance

Excellent temperature control (NG020A)

ORDERING INFORMATION		PRICE	PROMO
NG002	Mini Centrifuge, 2200 rpm/480g. With standard microtube, slide and strip tube rotor	244€	198 €
NG040	Plate Centrifuge, 6000 rpm/2000g. For 96- & 384-w. plates	943€	868 €
VX2	Vortexer Mini, 4.000rpm w/ adjust. speed	241€	198 €
NG020A	FastGene Mini Dry Bath Advance	313€	295 €

BIOTOP JL
www.biotop.fi

Itäinen Pitkäkatu 4 A
20520 TURKU
puh. (02) 241 0099
info@biotop.fi

Orders, technical support:

Matti Raitio
Risto Ranta
City-Lab Oy

puh. 050-361 1848 matti.rautio@biotop.fi
puh. 050-5280789 risto.ranta@biotop.fi
www.city-lab.fi

18.10.2022